- For post-curing of the parts, please select the parameter settings for your process. You can find more qualified printers on our website www.dentamid.dreve.de.

Post curing FotoDent® model2 beige-opaque	
PCU LED N_2 Rapid Shape DII / Asiga MAX UV Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	8 minutes, 80 % light power (vacuum)
FotoDent® flash Rapid Shape DII / Asiga MAX UV	8 minutes + 1 minute gas pre-flow / 12 minutes + 1 minute gas pre-flow
PCU 90	8 minutes (incl. 4 minutes gas pre-flow)
RS Cure	According to respective parameter set for FotoDent® model2
Post curing FotoDent® model2 black	
PCU LED N_2 Rapid Shape DII / Asiga MAX UV Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	12 minutes, 80 % light power (vacuum)
RS Cure	According to respective parameter set for FotoDent® model2
Post curing FotoDent® model2 white	
PCU LED N_2 Rapid Shape DII / Asiga MAX UV Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	6 minutes, 40 % light power (vacuum)
RS Cure	According to respective parameter set for FotoDent® model2

- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® model2. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hour with FotoDent® model2.

To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and/or colour variations of the FotoDent® model2 material.

Hazards identification

Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® model2 packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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FotoDent® model2 (EN)

Product description

Light-curing material for the fabrication of dental working models.

Processing

- During the application of FotoDent® model2 we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- FotoDent® model2 should be well shaken approx. 1 hour before use.
- Carefully pour FotoDent® model2 into the pre-determined container of the production unit.
- Please see further steps within printer manual.
- Remove all bubbles with a cleaned object.

Processing temperature	
Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	23 ± 3 °C
Rapid Shape DII printer	23 ± 3 °C
Asiga MAX UV	35 ± 3 °C

- Select the product specific parameter settings for FotoDent® model2 for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2 x 6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended. Alternatively, cleaning with isopropanol can be carried out in a separate beaker under pouring conditions for approx. 2 x 6 minutes. Recommendation: Beaker and magnetic stirrer, laboratory shaker or ultrasonic cleaner. In case of insufficient cleaning (such as shiny areas on the object after drying) repeat cleaning process.

Dreve Dentamid GmbH is not liable for any damages caused by improper application of the material.

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Dreve Dentamid GmbH haftet nicht für Schäden, die durch fehlerhafte Anwendung des Materials hervorgerufen werden.